

CITE-seq & Cell Hashing Protocol

For experiments involving Cell Hashing, use [cost per cell calculator](#) to plan experiments, determine number of hashes, number of cells to load, expected doublet rates (detected and undetected) and cost considerations:

Obtain cells

- Carefully count all cells to ensure accurate quantitation.
 - Make note of cell viability (>95%) and also include dead cells in the total cell count!
 - If you observe many dead cells, live cell enrichment is recommended!
 - Keep cell suspensions on ice (unless otherwise stated) at all times.

Multiplexing different samples

- Resuspend ~1-2 million cells of each sample in 100 µl Staining buffer (2%BSA/0.01%Tween, PBS) in 1.5 ml low bind tubes.
- Add 10 µl Fc Blocking reagent (FcX, BioLegend).
- Incubate for 10 minutes at 4°C.
- While cells are incubating in Fc Block, prepare antibody-pool using ~0.5 - 1 µg (or titrated amounts) of each **CITE-seq antibody**.
- Add **CITE-seq antibody pool** to each tube.
- Add 0.5ug of a unique **Cell Hashing antibody** to each tube.
- Incubate for 30 minutes at 4°C.
- Wash cells 3 times with 1 ml Staining buffer (2%BSA/0.01%Tween, PBS), spin at 4°C for 5 minutes at 350g.
- **Keep samples separated until shortly before the run.**
- Merge all samples at desired proportions in 1 ml Staining buffer, spin 5 minutes 350g at 4°C.
- Resuspend cells in PBS at appropriate concentration for downstream application.
(e.g. for **10x** ~500 cells/µl; for **Drop-seq** [~ 200 cells/µl]; for **super-loading** ~ 1,500 cells/µl or higher).
- Filter cells through 40 µm strainers (e.g. Flowmi cell strainer).
- **Verify cell concentration, single cellularity and high cell viability on hemocytometer after filtration.**
- **Immediately proceed to single cell run.**

Super-loading of one sample

- Resuspend ~1-2 million cells in 100 µl Staining buffer (2%BSA/0.01%Tween, PBS) in 1.5 ml low bind tube.
- Add 10 µl Fc Blocking reagent (FcX, BioLegend).
- Incubate for 10 minutes at 4°C.
- While cells are incubating in Fc Block, prepare antibody-pool using ~0.5 - 1 µg (or titrated amounts) of each **CITE-seq antibody**.
- Add **CITE-seq antibody pool** to cells.
- Incubate for 30 minutes at 4°C.
- Divide sample into 10 equal proportions into 1.5 ml low bind tubes (when using 10 Hashtags).
- Adjust volume to 50 µl with Staining buffer (2%BSA/0.01%Tween, PBS)
- Add 0.5ug of unique **Cell Hashing antibody** to each tube.
- Incubate for 20 minutes at 4°C.
- Wash cells 3 times with 1 ml Staining buffer (2%BSA/0.01%Tween, PBS), spin at 4°C for 5 minutes at 350g.
- **Keep samples separated until shortly before the run.**
- Merge cells from all tubes at equal ratios in 1 ml Staining buffer, spin 5 minutes 350g at 4°C.
- Resuspend cells in PBS at appropriate concentration for desired super-loading yield.
(~ 1,500 cells/µl or higher, use [cost per cell calculator](#)).
- Filter cells through 40 µm strainers (e.g. Flowmi cell strainer).
- **Verify cell concentration, single cellularity and high cell viability on hemocytometer after filtration.**
- **Immediately proceed to single cell run.**

Run **Drop-seq (Macosko et al., 2015)** or **10x Genomics single cell 3' (v2 or v3)** assay as described until before cDNA amplification.

When super-loading we recommend adding 0.5U/µl RNAse inhibitor to the single cell RT master mix.

At cDNA amplification step:**Add “additive” primers to cDNA PCR to increase yield of ADT and/or HTO products:****ADT PCR additive primer (0.2 μ M): 1 μ l (for 10x Genomics) or 0.4 μ l (for Drop-seq)****HTO PCR additive primer (0.1 μ M): 1 μ l (for 10x Genomics) or 0.4 μ l (for Drop-seq)**

Subtract the total volume of additive primer from the water added to the PCR reaction.

After cDNA amplification: Separate ADT / HTO-derived cDNAs (<180bp) and mRNA-derived cDNAs (>300bp).

- Perform SPRI selection to separate mRNA-derived and antibody-oligo-derived cDNAs.
- **DO NOT DISCARD SUPERNATANT FROM 0.6X SPRI. THIS CONTAINS THE ADTs and hashtags!**
 - Add 0.6X SPRI to cDNA reaction as described in 10x Genomics or Drop-seq protocol.
 - Incubate 5 minutes and place on magnet.
 - **Supernatant contains ADTs and Cell Hashtags.**
 - Beads contain full length mRNA-derived cDNAs.
- **mRNA-derived cDNA >300bp (beads fraction).**
 - Proceed with standard 10x or Drop-seq protocol for cDNA sequencing library preparation.
- **ADTs and Hashtags <180bp (supernatant fraction).**
 - Purify ADTs using two 2X SPRI purifications:
 - Add 1.4X SPRI to supernatant to obtain a final SPRI volume of 2X SPRI.
For example: For 10x Genomics, the cDNA PCR reaction volume is 100 μ l, 60 μ L (0.6X SPRI was added for capturing full length cDNAs, therefore add another 140 μ l of SPRI to the supernatant. The final volume will be 300 μ l containing 2X SPRI.
 - Transfer entire volume into a low-bind 1.5mL tube.
 - Incubate 10 minutes at room temperature.
 - Place tube on magnet and wait ~2 minutes until solution is clear.
 - Carefully remove and discard the supernatant.
 - Add 400 μ l 80% Ethanol to the tube without disturbing the pellet and stand for 30 seconds (only one Ethanol wash).
 - Carefully remove and discard the ethanol wash.
 - Centrifuge tube briefly and return it to magnet.
 - Remove and discard any remaining ethanol.
 - Resuspend in beads in 50 μ l water.
 - Perform another round of 2X SPRI purification by adding 100 μ l SPRI reagent directly onto resuspended beads.
 - Mix by pipetting and incubate 10 minutes at room temperature.
 - Place tube on magnet and wait ~2 minutes until solution is clear.
 - Carefully remove and discard the supernatant.
 - Add 200 μ l 80% Ethanol to the tube without disturbing the pellet and stand for 30 seconds (1st Ethanol wash).
 - Carefully remove and discard the ethanol wash.
 - Add 200 μ l 80% Ethanol to the tube without disturbing the pellet and stand for 30 seconds (2nd Ethanol wash).
 - Carefully remove and discard the ethanol wash.
 - Centrifuge tube briefly and return it to magnet.
 - Remove and discard any remaining ethanol and allow the beads to air dry for 2 minutes (do not over dry beads).
 - Resuspend beads in 90 μ l water.
 - Pipette mix vigorously and incubate at room temperature for 5 minutes.
 - Place tube on magnet and transfer clear supernatant into two PCR tubes.

- Amplify CITE-seq ADT sequencing library:
 - Prepare 100uL PCR reaction with purified ADTs:
 - 45 µl purified ADT/Hashtag fraction
 - 50 µl 2x KAPA Hifi PCR Master Mix.
 - 2.5 µl TruSeq Small RNA RPIx primer (containing i7 index) 10 µM.
 - 2.5 µl P5 oligo at 10 µM depending on application:
 - For Drop-seq use P5-SMART-PCR hybrid oligo.
 - For 10x use SI PCR oligo.
 - Cycling conditions:

95°C	3 min		
95°C	20 sec		
60°C	30 sec		~ 6-10 cycles
72°C	20 sec		
72°C	5 min		
- Amplify Cell Hashing HTO sequencing library:
 - Prepare 100uL PCR reaction with purified small fraction:
 - 45 µl purified ADT/Hashtag fraction
 - 50 µl 2x KAPA Hifi PCR Master Mix.
 - 2.5 µl TruSeq DNA D7xx_s primer (containing i7 index) 10 µM.
 - 2.5 µl P5 oligo at 10 µM depending on application:
 - For Drop-seq use P5-SMART-PCR hybrid oligo.
 - For 10x use SI PCR oligo.
 - Cycling conditions:

95°C	3 min		
95°C	20 sec		
64°C	30 sec		~ 8-12 cycles
72°C	20 sec		
72°C	5 min		
- Purify PCR products using 1.6X SPRI purification by adding 160 µl SPRI reagent.
 - Incubate 5 minutes at room temperature.
 - Place tube on magnet and wait 1 minute until solution is clear.
 - Carefully remove and discard the supernatant.
 - Add 200 µl 80% Ethanol to the tube without disturbing the pellet and stand for 30 seconds (first Ethanol wash).
 - Carefully remove and discard the ethanol wash.
 - Add 200 µl 80% Ethanol to the tube without disturbing the pellet and stand for 30 seconds (second Ethanol wash).
 - Carefully remove and discard the ethanol wash.
 - Centrifuge tube briefly and return it to magnet.
 - Remove and discard any remaining ethanol and allow the beads to air dry for 2 minutes.
 - Resuspend beads in 30 µl water.
 - Pipette mix vigorously and incubate at room temperature for 5 minutes.
 - Place tube on magnet and transfer clear supernatant to low-bind tube.
- ADT and Hashtag libraries are now ready to be sequenced.
 - Quantify libraries by standard methods (QuBit, BioAnalyzer, qPCR).
 - ADT and Hashtag libraries will be around 180 bp (Figure 1 and 2).

Oligonucleotide sequences:

CITE-seq antibody-oligos (ADTs):

CITE-seq antibody-oligos contain standard small TruSeq RNA read 2 sequences and can be amplified using Illumina's TruSeq Small RNA primer sets (RPIx – primers, see example RPI1 below). See example below with a 12nt barcode:

5' CCTTGGCACCCGAGAATTCCA **CATGATTGGCTC** BAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA***A*****A**

Hashtag barcoding antibody-oligos (HTOs):

Cell Hashing antibody-oligos contain standard TruSeq DNA read 2 sequences and can be amplified using truncated versions of Illumina's TruSeq DNA primer sets (see example D701_s below). See example below with a 12nt barcode:

5' GTGACTGGAGTTCAGACGCTGTGCTCTTCCGATCT **AGGACCATCCA** BAAAAAAAAAAAAAAAAAAAAAAAAAAAAAA***A*****A**

Oligos required for ADT and HTO library amplification:

- Drop-seq P5-SMART-PCR hybrid primer (for Drop-seq only)
5' AATGATACGGCGACCACCGAGATCTACACGCTGTCCGCGGAAGCAGTGGTATCAACGCAGAGT***A*****C**
 - 10x Genomics SI-PCR primer (for 10x Genomics Single Cell 3P (v2 and v3) and 5P kits)
5' AATGATACGGCGACCACCGAGATCTACACTCTTCCCTACACGACGCTC
 - ADT CITE-seq cDNA PCR additive primer
5' CCTTGGCACCCGAGAATTCC
 - HTO Cell Hashing cDNA PCR additive primer
5' GTGACTGGAGTTCAGACGCTGTGCTC
 - Illumina Small RNA RPI1 primer (for ADT amplification; i7 index 1, Oligonucleotide sequences © 2015 Illumina, Inc)
5' CAAGCAGAAGACGGCATACGAGAT**CTCGT**GATGTGACTGGAGTTCCTTGGCACCCGAGAATTCCA
 - Illumina TruSeq D701_s primer (for HTO amplification; i7 index 1, shorter than the original D701 Illumina sequence)
5' CAAGCAGAAGACGGCATACGAGAT**CGAGTAA**TGTGACTGGAGTTCAGACGCTGTC
- * Phosphorothioate bond
B C or G or T; not A nucleotide

Materials and Kits needed:

- Antibody-oligo conjugates for CITE-seq and/or Cell Hashing (e.g. [BioLegend Totalseq A](#))
- FC blocking reagent (e.g. BioLegend FcX)
- 8-strip PCR tubes, **emulsion safe (!)** (e.g. TempAssure PCR 8-strips, USA Scientific)
- Bioanalyzer chips and reagents (DNA High Sensitivity kit, Agilent)
- SPRiSelect reagent (GE Healthcare, B23317)
- Low-bind 1.5 mL tubes (Common lab suppliers)
- PCR Thermocycler (e.g. Bio-Rad, T100)
- Magnetic tube rack (e.g. Invitrogen, USA)
- Qubit (Invitrogen, USA)
- Hemocytometer (e.g. Fuchs Rosenthal)
- DMSO (Common lab suppliers)
- PBS (Common lab suppliers)
- Tween20 (Common lab suppliers)
- TE pH 8.0 (Common lab suppliers)
- BSA (DNase, RNase and protease free, e.g. VWR #0332-25G)
- Dead Cell Removal Kit (optional, e.g. Miltenyi)