

## ASAP-seq protocol

### 1. Cell staining

1. Obtain all single cell suspensions (filter if needed) and measure viability and density (if viability is <90% proceed with live cell enrichment and/or use best judgement depending on sample source / importance / rarity).
2. Resuspend ~1-2 million cells in 100  $\mu$ l Staining buffer (2%BSA, 0.01%Tween in PBS).
3. Add 10  $\mu$ l Fc Blocking reagent (FcX, BioLegend).
4. Incubate for 10 minutes at 4° C.
5. While cells are incubating in Fc Block, prepare antibody pool (panel or titrated amounts).
6. Add antibody-oligo pool to cells.
7. Incubate for 30 minutes at 4° C.
8. Wash cells 3 times with 1 mL Staining buffer spin 5 minutes 300g at 4° C.
9. Resuspend cells in 450  $\mu$ l PBS .

### 2. Cell fixation/lysis

1. Add 30  $\mu$ l 16% formaldehyde (1% f.c) and incubate 10min at room temperature, swirl, invert occasionally
2. Quench by adding glycine to 0.125M f.c.
3. Wash with 1x ice-cold PBS by filling up the tube, invert 5 times
4. Spin 5 minutes 400g at 4° C.
5. Discard supernatant and repeat wash with 1ml 1x ice-cold PBS
6. Spin 5 minutes 400g at 4° C, discard supernatant.
7. Resuspended cell pellet in 100  $\mu$ l chilled lysis buffer, mix by pipetting.
8. Incubate on ice for 3min (primary cells), 5min (cell lines)
9. Add 1 ml chilled wash buffer to the lysed cells, mix by pipetting
10. Spin 5 minutes 500g at 4° C.

(Only if intracellular staining is desired, otherwise skip to step 15):

11. Resuspend cell pellet in 50ul block buffer, ie 40ul intracellular wash buffer + 5ul FcX + 5ul monoblock, 15min on ice
12. Add 50ul of intracellular wash bf, containing indicated amount of conjugated intracellular markers, ice 30min
13. Wash with intracellular wash bf, 3x
14. Spin 5 minutes 500g at 4° C.

15. Remove supernatant, resuspend in 150  $\mu$ l 1x nuclei buffer (10x Genomics)
16. Filter through 40  $\mu$ m strainers (e.g. Flowmi cell strainer)
17. Count cells and adjust density according to 10x loading instructions.

### 3. Transposition and barcoding

18. Proceed according to 10x with the below modifications:

- 1) During the barcoding reaction (step 2.1) spike in 0.5  $\mu$ l of 1  $\mu$ M bridge oligo (there is no dead volume in the reaction, so final volume will be 65.5  $\mu$ l for v1 and 60.5  $\mu$ l for v1.1)
- 2) During GEM incubation (step 2.5), include a 5min incubation at 40°C at the beginning of the protocol, ie:

40°C 5 min	
72°C 5 min	
98°C 30sec	
98°C 10 sec	
59°C 30sec	repeat 11x (total 12 cycles)
72°C 1 min	
15°C hold	

(this extra step is not essential when using TSA products, but increases efficiency in TSB and especially TSC tag capture)

- 3) During silane bead elution (step 3.1o) add 43.5  $\mu$ l of Elution Solution I and subsequently recover ~43  $\mu$ l. Keep 3ul  $\mu$ l aside to use as input\* in the tag library PCR and with the remaining 40  $\mu$ l proceed to SPRI clean up as per protocol.
- 4) During SPRI cleanup (step 3.2d), save the supernatant. For the bead bound fraction proceed as per protocol. For the sup fraction, add 32  $\mu$ l SPRI, let bind for 5min. Collect beads on magnet, wash 2x with 80% EtOH, remove remaining ethanol and elute beads in 42  $\mu$ l EB (or more if multiple indexing reactions need to take place, ie TSA and TSB products have been used)\*. This can be combined with the 3  $\mu$ l left aside after the silane purification, as input in the ADT/HTO indexing reaction:

PCR:

50 $\mu$ l	2x KAPA mix
2.5 $\mu$ l	P5 10 $\mu$ M
2.5 $\mu$ l	RPxx (TSA family) or D7xx (TSB family) 10 $\mu$ M
3-45 $\mu$ l input fragments	
Water to 100 $\mu$ l total	

Program:

95°C	3 min	
95°C	20 sec	
60°C	30 sec	repeat ~13-15x (total 14-16 cycles)
72°C	20 sec	
72°C	5 min	

4°C hold

Expected product: 190 bp

\*note: you can use either as input in the tag indexing reaction. I find useful to combine when working with large antibody panels to increase input complexity. For simple hashing or small antibody panels you can rely on either alone.

**Buffers:**

<b>LYSIS</b>	LLL	OMNI	Stock
Tris-HCl (pH 7.5)	10mM	10mM	1M
NaCl	10mM	10mM	5M
MgCl <sub>2</sub>	3mM	3mM	1M
NP40 (IGEPAL)	0.1%	0.1%	10%
Digitonin	-	0.01%	5%
Tween	-	0.1%	10%
BSA	1%	1%	10%

<b>WASH</b>	Wash buffer
Tris-HCl (pH 7.5) 1M	10mM
NaCl 5M	10mM
MgCl <sub>2</sub> 1M	3mM
BSA 10%	1%

Intracellular staining buffer: Available from BioLegend. Part number 900002577

**Oligos:****BOA:**

TCGTCGGCAGCGTCAGATGTGTATAAGAGACAGNNNNNNNNNVTTTTTTTTTTTTTTTTTTTTTTTTTTTTTTTTTTTTTTT  
TT/3InvdT/

**BOB:**

TCGTCGGCAGCGTCAGATGTGTATAAGAGACAGTTGCTAGGACCGCCTTAAAGC/3InvdT/

**Example of an RPxx (TruSeq Small RNA handle, present in TSA tags) :**

CAAGCAGAAGACGGCATAACGAGATxxxxxxxGTGACTGGAGTTCCTTGGCACCCGAGAATTCCA

**Example of an D7xx (TruSeq DNA handle, present in TSB tags or TSA hashing):**

CAAGCAGAAGACGGCATAACGAGATxxxxxxxGTGACTGGAGTTCAGACGTGTGC

Replace “xxxxxxx” nucleotides with desired sample index